

JUNE
2023

siiora

Preserving Nature's Most Wonderful Creation



SCRUCAN[®]

Cannulated Compression System

Product Catalogue 2023



www.siiora.com



As an India Company we feel proud to showcase our capabilities in delivering excellence every time. With an international standard we continue to evolve our technology with pure Indian Values. We serve individual with our quality products. We believe in a non-ending relation with our associates.

Anuj Dureja
Director
Siora Surgicals Pvt. Ltd.

Founded in 1987, Siora has steadily grown to become a world leader of innovation in the manufacturing of orthopedic implants and equipment. Our manufacturing unit is set up according to WHO-GMP norms and quality system is audited to conform to ISO 13485:2016 and ISO 9001:2015 standard.

Great Quality

Every Siora products comes with an assurance of high quality and zero compromise in its production, performance, and safety.

Business Vision


To become the largest manufacturer of orthopaedic implants in the world, establishing the relation of the brand with excellence.

Guiding Principles

Our perseverance to achieve the highest standards have given us the leader status in all our businesses.



SCRUCAN[®] Cannulated Compression System




Cannulation	1.0mm
Hex Size	2.0mm
Corrosponding Screw Driver	2.0mm Tip
Corrosponding Guide Wire	0.9mm
Corrosponding cannulated Drill Bit	2.2mm


Cat No.	Length (mm)	Front Thread Length (A)mm	Shaft Length mm	Back Thread Length (B)mm
Titanium				
Ti.133.214	14	5.25	4.90	3.85
Ti.133.216	16	6.00	5.60	4.40
Ti.133.218	18	6.75	6.30	4.95
Ti.133.220	20	7.50	7.00	5.50
Ti.133.222	22	8.25	7.70	6.05
Ti.133.224	24	9.00	8.40	6.60
Ti.133.226	26	9.75	9.10	7.15
Ti.133.228	28	10.50	9.80	7.70
Ti.133.230	30	11.25	10.50	8.25
Ti.133.235	35	13.00	12.40	9.60
Ti.133.240	40	15.00	14.00	11.00
Ti.133.245	45	16.90	15.60	12.50
Ti.133.250	50	18.80	17.40	13.80

2.7mm/ 3.5mm Cannulated Screw - Titanium

Cat No.	Length (mm)
Ti.133.214	14
Ti.133.216	16
Ti.133.218	18
Ti.133.220	20
Ti.133.222	22
Ti.133.224	24
Ti.133.226	26
Ti.133.228	28
Ti.133.230	30
Ti.133.235	35
Ti.133.240	40
Ti.133.245	45
Ti.133.250	50



SCRUCAN[®] Cannulated Compression System



Cannulation	1.25mm
Hex Size	2.5mm
Corrosponding Screw Driver	2.5mm Tip
Corrosponding Guide Wire	0.9mm
Corrosponding cannulated Drill Bit	2.5mm

Cat No.	Length (mm)	Front Thread Length (A)mm	Shaft Length mm	Back Thread Length (B)mm
Titanium				
Ti.135.314	14	4.42	6.32	3.26
Ti.135.316	16	5.06	8.83	3.73
Ti.135.318	18	5.69	9.88	4.19
Ti.135.320	20	6.33	9.01	4.66
Ti.135.322	22	6.96	9.91	5.13
Ti.135.324	24	7.59	10.82	5.59
Ti.135.326	26	8.23	14.29	6.06
Ti.135.328	28	8.86	12.61	6.53
Ti.135.330	30	9.50	13.50	7.00
Ti.135.335	35	11.08	19.24	8.16
Ti.135.340	40	12.66	18.01	9.33
Ti.135.345	45	14.24	20.27	10.49
Ti.135.350	50	15.83	22.51	11.66

3.5mm / 4.5mm Cannulated Screw - Titanium

Cat No.	Length (mm)
Ti.135.314	14
Ti.135.316	16
Ti.135.318	18
Ti.135.320	20
Ti.135.322	22
Ti.135.324	24
Ti.135.326	26
Ti.135.328	28
Ti.135.330	30
Ti.135.335	35
Ti.135.340	40
Ti.135.345	45
Ti.135.350	50



Instruments for **SCRUCAN®** Cannulated Compression System

Graphic Set for **SCRUCAN®**
Cannulated Compression System

FOR 2.7/3.5 & 3.5/4.5MM



Cat No. IS.133	Description	Units
TI.133.214	2.7mm/3.5mm Cannulated Screw - Titanium Length 14mm	2
TI.133.216	2.7mm/3.5mm Cannulated Screw - Titanium Length 16mm	2
TI.133.218	2.7mm/3.5mm Cannulated Screw - Titanium Length 18mm	2
TI.133.220	2.7mm/3.5mm Cannulated Screw - Titanium Length 20mm	2
TI.133.222	2.7mm/3.5mm Cannulated Screw - Titanium Length 22mm	2
TI.133.224	2.7mm/3.5mm Cannulated Screw - Titanium Length 24mm	2
TI.133.226	2.7mm/3.5mm Cannulated Screw - Titanium Length 26mm	2
TI.133.228	2.7mm/3.5mm Cannulated Screw - Titanium Length 28mm	2
TI.133.230	2.7mm/3.5mm Cannulated Screw - Titanium Length 30mm	2
TI.133.235	2.7mm/3.5mm Cannulated Screw - Titanium Length 35mm	2
TI.133.240	2.7mm/3.5mm Cannulated Screw - Titanium Length 40mm	2
TI.133.245	2.7mm/3.5mm Cannulated Screw - Titanium Length 45mm	2
TI.133.250	2.7mm/3.5mm Cannulated Screw - Titanium Length 50mm	2

Continues...

Instruments[®] for

SCRUCAN[®]

Cannulated Compression System


Graphic Set for **SCRUCAN[®]** Cannulated Compression System

FOR 2.7/3.5 & 3.5/4.5MM

...Continuation of page 5

Cat No. IS.133	Description	Units
TI.135.314	3.5mm/4.5mm Cannulated Screw - Titanium Length 14mm	2
TI.135.316	3.5mm/4.5mm Cannulated Screw - Titanium Length 16mm	2
TI.135.318	3.5mm/4.5mm Cannulated Screw - Titanium Length 18mm	2
TI.135.320	3.5mm/4.5mm Cannulated Screw - Titanium Length 20mm	2
TI.135.322	3.5mm/4.5mm Cannulated Screw - Titanium Length 22mm	2
TI.135.324	3.5mm/4.5mm Cannulated Screw - Titanium Length 24mm	2
TI.135.326	3.5mm/4.5mm Cannulated Screw - Titanium Length 26mm	2
TI.135.328	3.5mm/4.5mm Cannulated Screw - Titanium Length 28mm	2
TI.135.330	3.5mm/4.5mm Cannulated Screw - Titanium Length 30mm	2
TI.135.335	3.5mm/4.5mm Cannulated Screw - Titanium Length 35mm	2
TI.135.340	3.5mm/4.5mm Cannulated Screw - Titanium Length 40mm	2
TI.135.345	3.5mm/4.5mm Cannulated Screw - Titanium Length 45mm	2
TI.135.350	3.5mm/4.5mm Cannulated Screw - Titanium Length 50mm	2
802.114	Direct Measuring Device for 2.7mm/3.5mm/4.5mm Cannulated Screws	1
704.030	Cannulated Screw Driver 2.0mm Tip - Silicon Handle	1
704.012	Cannulated Screw Driver 2.5mm Tip - Silicon Handle	1
XP.802.009	Cannulated drill Bit Plain Shank - Dia. 2.2mm X 100 mm Length (for 2.7/3.5mm Cannulated Screws)	1
XP.802.014	Cannulated drill Bit Plain Shank - Dia. 2.5mm X 100mm Length (for 3.5/4.5mm Cannulated Screws)	1
802.121	Threaded Guide Wire 0.9mm X 150mm Length	5
802.102	Guide Wire Sleeve for 0.9mm Guide Wire - Outer Dia. 2.2mm	1
802.116	Guide Wire Sleeve for 0.9mm Guide Wire - Outer Dia. 2.5mm	1
802.104	Drill Sleeve for 2.2mm Drill Bit	1
802.106	Drill Sleeve for 2.5mm Drill Bit	1
802.107	Handle for 0.9 Guide Wire, 2.2mm & 2.5mm Drill Bit Sleeves	1
401.051	Sterilization Tube 160mm Length for Guide Wire 0.9mm	1
400.132	Graphic Aluminium Box of Implants & Instruments For 2.7/3.5 & 3.5/4.5mm Cannulated Screws	1

Instruments for **SCRUCAN[®]**



802.114

Direct Measuring Device for
2.7/3.5/4.5mm Cannulated Screws



704.030

Cannulated Screw Driver
2.0mm Tip - Silicon Handle



704.012

Cannulated Hexagonal Screw Driver
2.5mm Tip - Silicon Handle




XP802.009

Cannulated drill Bit Plain Shank - Dia.
2.2mm X 100 mm Length
(for 2.7/3.5mm Cannulated Screws)



XP802.014


Cannulated drill Bit Plain Shank - Dia.
2.5mm X 100mm Length
(for 3.5/4.5mm Cannulated Screws)



802.121


Threaded Guide Wire
0.9mm X 150mm Length

Instruments for **SCRUCAN®**



802.102

Guide Wire Sleeve for
0.9mm Guide Wire -
Outer Dia. 2.2mm



802.116

Guide Wire Sleeve for
0.9mm Guide Wire -
Outer Dia. 2.5mm




802.106

Drill Sleeve for 2.5mm Drill Bit




802.104

Drill Sleeve for 2.2mm Drill Bit



802.107


Handle for 0.9 Guide Wire,
2.2mm & 2.5mm Drill Bit Sleeve



401.051

Sterilization Tube 160mm Length
For Guide Wire 0.9mm

SCRUCAN[®] Cannulated Compression System



Corrosponding Screw Driver	2.5mm Tip
Corrosponding Guide Wire	1.2mm
Corrosponding Cannulated Drill Bit	3.5mm
Cannulation	1.7mm
Hex Size	2.5mm


Cat No.	Length (mm)	Front Thread Length (A)mm	Shaft Length mm	Back Thread Length (B)mm
Titanium				
Ti.134.230	30	7.98	16.36	5.66
Ti.134.235	35	9.31	19.08	6.61
Ti.134.240	40	10.64	21.81	7.55
Ti.134.245	45	12.0	24.5	8.5
Ti.134.250	50	13.3	27.26	9.44
Ti.134.255	55	14.63	29.99	10.38
Ti.134.260	60	15.96	32.72	11.32
Ti.134.265	65	17.29	35.44	12.27
Ti.134.270	70	18.62	38.17	13.21
Ti.134.275	75	19.95	40.89	14.16
Ti.134.280	80	21.28	43.62	15.10

4.5mm/ 5.5mm Cannulated Screw - Titanium

Cat No.	Length (mm)
Ti.134.230	30
Ti.134.235	35
Ti.134.240	40
Ti.134.245	45
Ti.134.250	50
Ti.134.255	55
Ti.134.260	60
Ti.134.265	65
Ti.134.270	70
Ti.134.275	75
Ti.134.280	80



SCRUCAN[®] Cannulated Compression System



Corrosponding Screw Driver	3.5mm Tip
Corrosponding Guide Wire	1.5mm
Corrosponding Cannulated Drill Bit	4.0mm
Cannulation	2.2mm
Hex Size	3.5mm

5.5mm/ 6.5mm Cannulated Screw - Titanium

Cat No.	Length (mm)
TI.136.250	50
TI.136.255	55
TI.136.260	60
TI.136.265	65
TI.136.270	70
TI.136.275	75
TI.136.280	80
TI.136.285	85
TI.136.290	90
TI.136.295	95
TI.136.300	100

Cat No.	Length (mm)	Front Thread Length (A)mm	Shaft Length mm	Back Thread Length (B)mm
Titanium				
TI.136.250	50	12.9	28.2	8.9
TI.136.255	55	14.1	31.0	9.8
TI.136.260	60	15.4	33.9	10.7
TI.136.265	65	16.7	36.7	11.6
TI.136.270	70	18.0	39.5	12.5
TI.136.275	75	19.3	42.3	13.4
TI.136.280	80	20.6	45.1	14.3
TI.136.285	85	21.9	48.0	15.2
TI.136.290	90	23.1	50.8	16.1
TI.136.295	95	24.4	53.6	17.0
TI.136.300	100	25.7	56.4	17.9




Instruments for **SCRUCAN®** Cannulated Compression System

Graphic Set for **SCRUCAN®**
Cannulated Compression
System

FOR 4.5/5.5 & 5.5/6.5MM

★★



Cat No. IS.134	Description	Units
TI.134.230	4.5mm/5.5mm Cannulated Screw - Titanium Length 30mm	2
TI.134.235	4.5mm/5.5mm Cannulated Screw - Titanium Length 35mm	2
TI.134.240	4.5mm/5.5mm Cannulated Screw - Titanium Length 40mm	2
TI.134.245	4.5mm/5.5mm Cannulated Screw - Titanium Length 45mm	2
TI.134.250	4.5mm/5.5mm Cannulated Screw - Titanium Length 50mm	2
TI.134.255	4.5mm/5.5mm Cannulated Screw - Titanium Length 55mm	2
TI.134.260	4.5mm/5.5mm Cannulated Screw - Titanium Length 60mm	2
TI.134.265	4.5mm/5.5mm Cannulated Screw - Titanium Length 65mm	2
TI.134.270	4.5mm/5.5mm Cannulated Screw - Titanium Length 70mm	2
TI.134.275	4.5mm/5.5mm Cannulated Screw - Titanium Length 75mm	2
TI.134.280	4.5mm/5.5mm Cannulated Screw - Titanium Length 80mm	2

Continues...

Instruments ^{for} **SCRUCAN[®]**

Cannulated Compression System

Graphic Set for **SCRUCAN[®]** Cannulated Compression System

FOR 4.5/5.5 & 5.5/6.5MM

...Continuation of page 11

Cat No. IS.134	Description	Units
TI.136.250	5.5mm/6.5mm Cannulated Screw - Titanium Length 50mm	2
TI.136.255	5.5mm/6.5mm Cannulated Screw - Titanium Length 55mm	2
TI.136.260	5.5mm/6.5mm Cannulated Screw - Titanium Length 60mm	2
TI.136.265	5.5mm/6.5mm Cannulated Screw - Titanium Length 65mm	2
TI.136.270	5.5mm/6.5mm Cannulated Screw - Titanium Length 70mm	2
TI.136.275	5.5mm/6.5mm Cannulated Screw - Titanium Length 75mm	2
TI.136.280	5.5mm/6.5mm Cannulated Screw - Titanium Length 80mm	2
TI.136.285	5.5mm/6.5mm Cannulated Screw - Titanium Length 85mm	2
TI.136.290	5.5mm/6.5mm Cannulated Screw - Titanium Length 90mm	2
TI.136.295	5.5mm/6.5mm Cannulated Screw - Titanium Length 95mm	2
TI.136.300	5.5mm/6.5mm Cannulated Screw - Titanium Length 100mm	2
802.115	Direct Measuring Device for 4.5mm/5.5mm/6.5mm Cannulated Screws	1
802.122	Threaded Guide Wire 1.2mm, 200mm Length	5
704.012	Cannulated Hexagonal Screw Driver 2.5mm Tip - Silicon Handle	1
703.016	Cannulated Hexagonal Screw Driver 3.5mm Tip - Silicon Handle	1
XP.802.011	Cannulated Drill Bit Plain Shank - Dia. 3.5mm X 200mm Length	1
XP.802.012	Cannulated Drill Bit Plain Shank - Dia. 4.0mm X 200mm Length	1
	(for 4.5/5.5mm Cannulated Screws)	
	(for 5.5/6.5mm Cannulated Screws)	
802.013	Guide Wire 1.5mm X 200mm Length	5
802.108	Guide Wire Sleeve for 1.2mm & 1.5mm Guide Wire - Outer Dia. 3.5mm	1
802.117	Guide Wire Sleeve for 1.2mm & 1.5mm Guide Wire - Outer Dia. 4.0mm	1
802.110	Drill Sleeve 8.0mm for 3.5mm Drill Bit	1
802.112	Drill Sleeve 8.0mm for 4.0mm Drill Bit	1
802.113	Handle for 1.2mm & 1.5mm Guide Wire, 3.5mm & 4.0mm Drill Bit Sleeves	1
401.061	Sterilization Tube 210mm Length, for Guide Wire 1.5mm & 1.2mm	1
400.133	Graphic Aluminium Box of Implants & Instruments For 4.5/5.5 & 5.5/6.5mm Cannulated Screws	1

Instruments for **SCRUCAN[®]**




802.115

Direct Measuring Device for
4.5/5.5mm/6.5mm Cannulated Screws



703.016

Large Cannulated Hexagonal Screw Driver
3.5mm Tip - Silicon Handle



704.012

Cannulated Hexagonal Screw Driver
2.5mm Tip - Silicon Handle



XP802.011

Cannulated Drill Bit Plain Shank - Dia.
3.5mm X 200mm Length
(for 4.5/5.5mm Cannulated Screws)




XP802.012

Cannulated Drill Bit
Plain Shank - Dia.
4.0mm X 200mm Length
(for 5.5/6.5mm Cannulated Screws)



802.013

Guide Wire
1.5mm X 200mm Length



802.122

Threaded Guide Wire
1.2mm X 200mm Length

Instruments for **SCRUCAN®**



802.108

Guide Wire Sleeve for
1.2mm & 1.5mm
Guide Wire - Outer Dia. 3.5mm



401.061

Sterilization Tube 210mm Length
For Guide Wire 1.5mm & 1.2mm



802.117

Guide Wire Sleeve for 1.2mm & 1.5mm
Guide Wire - Outer Dia. 4.0mm



802.110

Drill Sleeve 8.0mm for 3.5mm Drill Bit



802.112

Drill Sleeve 8.0mm for 4.0mm Drill Bit



802.113

Handle for 1.2mm & 1.5mm Guide Wire,
3.5mm & 4.0mm Drill Bit Sleeve

Issue No. - 6

Revision No. - 6

UP DATED : 06/2023



Siora Surgicals Private Limited

Manufacturer of :

ORTHOPAEDIC IMPLANTS & INSTRUMENTS

Regd. & Corporate Office

WZ-1, Phool Bagh, Rampura, Delhi - 110 035, INDIA. Phone: +91 11 28316116, 17, 18

Manufacturing Unit:

1792, HSIIDC Indl. Estate, RAI District, Sonapat, Haryana - 131029, INDIA

Website:

www.siiora.com

DISCLAIMER: This catalog is intended exclusively for the reference of Physicians. Any information about the products and procedures listed here is general in nature and may not be construed as any medical advice, whatsoever. Information so provided herein therefore does not qualify to be any diagnostic and / or therapeutic opinion with regard to any individual medical case and a patient may be examined and advised accordingly. The representations thus made in the contents of the catalog are not an alternative to any such examination.