



LOCKING PLATE SYSTEM

PRODUCT CATALOG

3.5mm Pediatric Eight-plate



3.5mm Cannulated Screw Full-threaded



4.5 Pediatric Hinge Locking plate



4.5mm Pediatric Eight-plate



4.5mm Pediatric Plate



4.5mm Cannulated Screw Full-threaded



REF	Holes	Length
01.1261-01	2 holes*10mm	10mm
01.1261-02	2 holes*12mm	12mm
01.1261-03	2 holes*16mm	16mm

Use 3.5mm Cannulated Screw Full-threaded

REF	Specification	REF	Specification
01.1261-04	3.5*10mm	01.1261-12	3.5*26mm
01.1261-05	3.5*12mm	01.1261-13	3.5*28mm
01.1261-06	3.5*14mm	01.1261-14	3.5*30mm
01.1261-07	3.5*16mm	01.1261-15	3.5*32mm
01.1261-08	3.5*18mm	01.1261-16	3.5*35mm
01.1261-09	3.5*20mm	01.1261-17	3.5*40mm
01.1261-10	3.5*22mm	01.1261-18	3.5*45mm
01.1261-11	3.5*24mm	01.1261-19	3.5*50mm

REF	Holes	Length
01.1264-01	2 holes*12mm	12mm
01.1264-02	2 holes*16mm	16mm

Use 4.5mm Cannulated Screw Full-threaded

REF	Holes	Length
01.1262-01	2 holes*12mm	12mm
01.1262-02	2 holes*16mm	16mm
01.1262-03	2 holes*20mm	20mm
01.1262-04	2 holes*25mm	25mm

Use 4.5mm Cannulated Screw Full-threaded

REF	Holes	Length
01.1263-01	4 holes*12mm	12mm
01.1263-02	4 holes*16mm	16mm
01.1263-03	4 holes*20mm	20mm

Use 4.5mm Cannulated Screw Full-threaded

REF	Specification	REF	Specification
01.1262-05	4.5*10mm	01.1262-13	4.5*26mm
01.1262-06	4.5*12mm	01.1262-14	4.5*28mm
01.1262-07	4.5*14mm	01.1262-15	4.5*30mm
01.1262-08	4.5*16mm	01.1262-16	4.5*32mm
01.1262-09	4.5*18mm	01.1262-21	4.5*34mm
01.1262-10	4.5*20mm	01.1262-22	4.5*36mm
01.1262-11	4.5*22mm	01.1262-23	4.5*38mm
01.1262-12	4.5*24mm	01.1262-18	4.5*40mm

3.5mm Pediatric Sliding Barrel Locking Plate



3.5mm Pediatric Sliding Straight Locking Plate



2.7mm LCP Pediatric Hip Plate



3.5mm LCP Pediatric Hip Plate



2.7mm LCP Pediatric Straight Plate



REF	Holes	Length
01.1273-01	4 holes L	80mm
01.1273-02	4 holes R	80mm

Use 3.5mm Locking Screw

REF	Holes	Length
01.1274-01	5 holes	60mm
01.1274-02	7 holes	80mm
01.1274-03	9 holes	100mm
01.1274-04	11 holes	120mm

Use 3.5mm Locking Screw

REF	Holes	Length
01.7003-01	3 holes*100°	48.2mm
01.7003-02	3 holes*110°	48.2mm
01.7003-03	4 holes*120°	58.2mm
01.7003-04	6 holes*120°	78.2mm
01.7003-05	3 holes*130°	48.2mm
01.7003-06	10 holes*130°	118.2m

Use 2.7mm Locking Screw, 2.7mm Cortical Screw

REF	Holes	Length
01.7004-01	3 holes*90°	75mm
01.7004-02	4 holes*90°	87mm
01.7004-03	5 holes*90°	99mm
01.7004-04	6 holes*90°	111mm
01.7004-05	7 holes*90°	123mm
01.7004-06	8 holes*90°	135mm
01.7004-07	9 holes*90°	147mm
01.7004-08	3 holes*100°	75mm
01.7004-09	3 holes*110°	75mm
01.7004-10	3 holes*120°	62mm
01.7004-11	4 holes*120°	75mm
01.7004-12	6 holes*120°	101mm
01.7004-13	8 holes*120°	127mm
01.7004-14	10 holes*120°	153mm
01.7004-15	3 holes*130°	75mm
01.7004-16	4 holes*130°	87mm
01.7004-17	5 holes*130°	99mm
01.7004-18	7 holes*130°	123mm
01.7004-19	9 holes*130°	147mm
01.7004-20	3 holes*150°	64mm

Use 3.5mm Locking Screw, 3.5mm Cortical Screw

REF	Holes	Length
01.7010-01	5 holes	49.5mm
01.7010-02	6 holes	58.5mm
01.7010-03	7 holes	67.5mm
01.7010-04	8 holes	76.5mm
01.7010-05	9 holes	85.5mm
01.7010-06	10 holes	94.5mm

Use 2.7mm Locking Screw, 2.7mm Cortical Screw

5.0mm LCP Pediatric Hip Plate



REF	Holes	Length
01.7005-01	3 holes*90°	96mm
01.7005-02	4 holes*90°	112mm
01.7005-03	5 holes*90°	122mm
01.7005-04	7 holes*90°	220mm
01.7005-05	9 holes*90°	252mm
01.7005-06	3 holes*100°	101mm
01.7005-07	3 holes*110°	100mm
01.7005-08	3 holes*120°	99mm
01.7005-09	4 holes*120°	115mm
01.7005-10	6 holes*120°	147mm
01.7005-11	8 holes*120°	179mm
01.7005-12	10 holes*120°	211mm
01.7005-13	3 holes*130°	96mm
01.7005-14	4 holes*130°	112mm
01.7005-15	5 holes*130°	122mm

Use 5.0mm Locking Screw, 4.5mm Cortical Screw

2.7mm Cortex Screw, T8 Stardrive Self-tapping



REF	Specification	REF	Specification
030390010	2.7*10mm	030390022	2.7*22mm
030390012	2.7*12mm	030390024	2.7*24mm
030390014	2.7*14mm	030390026	2.7*26mm
030390016	2.7*16mm	030390028	2.7*28mm
030390018	2.7*18mm	030390030	2.7*30mm
030390020	2.7*20mm		

2.7mm Locking Screw, T8 Stardrive Self-tapping



REF	Specification	REF	Specification
030340010	2.7*10mm	030340026	2.7*26mm
030340012	2.7*12mm	030340028	2.7*28mm
030340014	2.7*14mm	030340030	2.7*30mm
030340016	2.7*16mm	030340032	2.7*32mm
030340018	2.7*18mm	030340034	2.7*34mm
030340020	2.7*20mm	030340036	2.7*36mm
030340022	2.7*22mm	030340038	2.7*38mm
030340024	2.7*24mm	030340040	2.7*40mm

3.5mm Locking Screw (Stardrive)



REF	Specification	REF	Specification
5100-3901	3.5*12mm	5100-3917	3.5*44mm
5100-3902	3.5*14mm	5100-3918	3.5*46mm
5100-3903	3.5*16mm	5100-3919	3.5*48mm
5100-3904	3.5*18mm	5100-3920	3.5*50mm
5100-3905	3.5*20mm	5100-3921	3.5*52mm
5100-3906	3.5*22mm	5100-3922	3.5*54mm
5100-3907	3.5*24mm	5100-3923	3.5*56mm
5100-3908	3.5*26mm	5100-3924	3.5*58mm
5100-3909	3.5*28mm	5100-3925	3.5*60mm
5100-3910	3.5*30mm	5100-3926	3.5*65mm
5100-3911	3.5*32mm	5100-3927	3.5*70mm
5100-3912	3.5*34mm	5100-3928	3.5*75mm
5100-3913	3.5*36mm	5100-3929	3.5*80mm
5100-3914	3.5*38mm	5100-3930	3.5*85mm
5100-3915	3.5*40mm	5100-3931	3.5*90mm
5100-3916	3.5*42mm	5100-3932	3.5*95mm

5.0mm Locking Screw (Stardrive)



REF	Specification	REF	Specification
5100-4001	5.0*22mm	5100-4015	5.0*50mm
5100-4002	5.0*24mm	5100-4016	5.0*52mm
5100-4003	5.0*26mm	5100-4017	5.0*54mm
5100-4004	5.0*28mm	5100-4018	5.0*56mm
5100-4005	5.0*30mm	5100-4019	5.0*58mm
5100-4006	5.0*32mm	5100-4020	5.0*60mm
5100-4007	5.0*34mm	5100-4021	5.0*65mm
5100-4008	5.0*36mm	5100-4022	5.0*70mm
5100-4009	5.0*38mm	5100-4023	5.0*75mm
5100-4010	5.0*40mm	5100-4024	5.0*80mm
5100-4011	5.0*42mm	5100-4025	5.0*85mm
5100-4012	5.0*44mm	5100-4026	5.0*90mm
5100-4013	5.0*46mm	5100-4027	5.0*95mm
5100-4014	5.0*48mm		

3.5mm Cortical Screw (Stardrive)





REF	Specification	REF	Specification
5100-4001	5.0*22mm	5100-4015	5.0*50mm
5100-4002	5.0*24mm	5100-4016	5.0*52mm
5100-4003	5.0*26mm	5100-4017	5.0*54mm
5100-4004	5.0*28mm	5100-4018	5.0*56mm
5100-4005	5.0*30mm	5100-4019	5.0*58mm
5100-4006	5.0*32mm	5100-4020	5.0*60mm
5100-4007	5.0*34mm	5100-4021	5.0*65mm
5100-4008	5.0*36mm	5100-4022	5.0*70mm
5100-4009	5.0*38mm	5100-4023	5.0*75mm
5100-4010	5.0*40mm	5100-4024	5.0*80mm
5100-4011	5.0*42mm	5100-4025	5.0*85mm
5100-4012	5.0*44mm	5100-4026	5.0*90mm
5100-4013	5.0*46mm	5100-4027	5.0*95mm
5100-4014	5.0*48mm		

4.5mm Cortical Screw (Stardrive)





REF	Specification	REF	Specification
5100-4201	4.5*22mm	5100-4211	4.5*42mm
5100-4202	4.5*24mm	5100-4212	4.5*44mm
5100-4203	4.5*26mm	5100-4213	4.5*46mm
5100-4204	4.5*28mm	5100-4214	4.5*48mm
5100-4205	4.5*30mm	5100-4215	4.5*50mm
5100-4206	4.5*32mm	5100-4216	4.5*52mm
5100-4207	4.5*34mm	5100-4217	4.5*54mm
5100-4208	4.5*36mm	5100-4218	4.5*56mm
5100-4209	4.5*38mm	5100-4219	4.5*58mm
5100-4210	4.5*40mm	5100-4220	4.5*60mm

Pediatric Hinge Locking Plate Instrument Set 3200-18



NO.	REF	Description	Qty.
1	3200-1801	Screwdriver SW2.5	1
2	3200-1802	Cannulated Screwdriver SW2.5	1
3	3200-1803	K-wire 1.5*150	1
4	3200-1804	K-wire 1.5*150	1
5	3200-1805	K-wire 1.5*150	1
6	3200-1806	K-wire 1.5*150	1
7	3200-1807	Target Mould Houlder	1
8	3200-1808	Target Mould Houlder	1
9	3200-1809	K-wire Depth Gauge 0-60mm	1
10	3200-1810	Screw Holding Forcep	1
11	3200-1811	Hign Plate Holding Forcep	1
1.2	3200-1812	Drill Sleeve 4.5mm	1
13	3200-1813	Drill Sleeve 4.5mm	1
14	3200-1814	Hex Wrench	1
15	3200-1815	Plate & Screw Box	1
16	3200-1816	Target Mould 2 holes*10mm	1
17	3200-1817	Target Mould 2 holes*12mm	1
18	3200-1818	Target Mould 2 holes*18mm	1
19	3200-1819	Target Mould 2 holes*20mm	1
20	3200-1820	Target Mould 2 holes*25mm	1
21	3200-1821	Target Mould 4 holes*16mm	1
22	3200-1822	Target Mould 4 holes*20mm	2
23	3200-1823	Target Mould 4 holes*12mm	1
24	3200-1824	1.5mm Locking Sleeve	1
25	3200-1825	1.5mm Locking Sleeve	1
26	3200-1826	2.5mm Drill Sleeve	1
27	3200-1827	2.5mm Drill Sleeve	1
28	3200-1828	Sleeve	1
29	3200-1829	Sleeve	1
30	3200-1830	10mm Target Mould Holder	1
31	3200-1831	10mm Target Mould Holder	1
32	3200-1832	1.5mm Wire Sleeve	1
33	3200-1833	1.5mm Wire Sleeve	1
34	3200-1834	3.2mm Drill Sleeve	1
35	3200-1835	1.2mm Wire Sleeve	1
36	3200-1836	1.2mm Wire Sleeve	1
37	3200-1837	3.2mm Drill Sleeve	1
38	3200-1838	Cannulated Drill Bit 3.2*200	1
39	3200-1839	Cannulated Drill Bit 3.2*200 Blcok	1
40	3200-1840	Cannulated Drill Bit 2.5*200	1
41	3200-1841	Cannulated Drill Bit 2.5*200 Blcok	1
42	3200-1842	K-wire 1.2*150	1
43	3200-1843	K-wire 1.2*150	1
44	3200-1844	Aluminium Box	1

2.7/3.5/5.0 Pediatric Osteotomy Locking Plate Instrument Set 3200-22



NO.	REF	Description	Qty.
1	3200-2201	4.2mm Locking Sleeve	2
2	3200-2202	4.2mm Locking Sleeve	1
3	3200-2203	3.0mm Locking Sleeve	2
4	3200-2204	3.0mm Locking Sleeve	1
5	3200-2205	Tap HA2.7*70	1
6	3200-2206	Tap HA4.5*150	1
7	3200-2207	Tap HA3.5*140	1
8	3200-2208	T15 Star Screwdriver	2
9	3200-2209	T20 Star Screwdriver	2
10	3200-2210	T8 Star Screwdriver	2
11	3200-2211	AO Drill Bit 3.2*150	2
12	3200-2212	AO Drill Bit 2.0*200	2
13	3200-2213	AO Drill Bit 2.5*200	2
14	3200-2214	K-wire Depth Gauge	1
15	3200-2215	K-wire 2.0*230	4
16	3200-2216	K-wire 2.5*250	4
17	3200-2217	Threaded K-wire 2.0*230	4
18	3200-2218	Threaded K-wire 2.5*300	4
19	3200-2219	AO Drill Bit 4.2*300	2
20	3200-2220	AO Drill Bit 2.9*250	2
21	3200-2221	2.7mm Locking Sleeve	2
22	3200-2222	2.7mm Locking Sleeve	1
23	3200-2223	2.7mm Sleeve	1
24	3200-2224	2.7mm Sleeve	1
25	3200-2225	2.7 Locating Device	2
26	3200-2226	3.5/5.0 Locating Device	2
27	3200-2227	3.5/5.0 Angle Device 80-160°	1
28	3200-2228	2.7 Angle Device 80-160°	1
29	3200-2229	Triangle Location Plate 90°/50°/40°	1
30	3200-2230	Triangle Location Plate 80°/70°/30°	1
31	3200-2231	Triangle Location Plate 100°/60°/20°	1
32	3200-2232	T-handle	1
33	4200-0119	Oblique Reduction Forcep (170mm)	1
34	4200-0210	Oblique Reduction Forcep (230mm)	1
35	4200-0117	Self-centering Bone Holding Forcep (190mm)	1
36	4200-0214	Self-centering Bone Holding Forcep (270mm)	1
37	3200-2233	Guider 2.0/2.7	1
38	3200-2234	Guider 2.5/3.5	1
39	3200-2235	Guider 3.2/4.5	1
40	3200-2236	Straight Handle	1
41	3200-2237	Torque Wrench 1.5 N.m	1
42	3200-2238	2.7 Osteotome Locator	1
43	3200-2239	3.5/5.0 Osteotome Locator	1
44	3200-2240	2.7 Target Device	1
45	3200-2241	3.5 Target Device	1
46	3200-2242	5.0 Target Device	1
47	3200-2243	2.7 Depth Gauge	1
48	3200-2244	3.5 Depth Gauge	1
49	3200-2245	5.0 Depth Gauge	1
50	3200-2246	Torque Wrench 4.0 N.m	1
51	3200-2247	Torque Wrench 0.8 N.m	1
52	3200-2248	Aluminium Box	1

1.5S Mini Locking Plate



REF	Holes	Length
021260004	4 holes	23mm
021260006	6 holes	36mm

Thickness: 1.0mm, Width: 4.5mm

1.5S Mini Reconstruction Locking Plate



REF	Holes	Length
02127004	4 holes	20mm
02127006	6 holes	30mm
02127008	8 holes	40mm
02127010	10 holes	50mm

Thickness: 1.0mm, Width: 4.5mm

1.5S Mini T Locking Plate (3 holes)



REF	Holes	Length
021280004	4 holes	25mm
021280006	6 holes	35mm
021280008	8 holes	45mm

Thickness: 1.0mm, Width: 4.5mm

1.5S Mini T Locking Plate (4 holes)



REF	Holes	Length
020980004	4 holes	25mm
020980006	6 holes	35mm
020980008	8 holes	45mm

Thickness: 1.0mm, Width: 4.5mm

1.5S Mini Y Locking Plate



REF	Holes	Length
020990004	4 holes	26mm
020990006	6 holes	36mm
020990008	8 holes	46mm

Thickness: 1.0mm, Width: 4.5mm

1.5S Mini Condylar Locking Plate



REF	Holes	Length
021000004	4 holes	26mm
021000006	6 holes	36mm
021000008	8 holes	46mm

Thickness: 1.0mm, Width: 4.5mm

1.5S Mini Strut Locking Plate
(Oblique-angled)

REF	Holes
021011000	8 holes L
021012000	8 holes R

Thickness: 1.0mm

2.0S Mini Locking Plate



REF	Holes	Length
021020004	4 holes	31mm
021020005	5 holes	38mm
021020006	6 holes	45mm
021020008	8 holes	59mm

Thickness: 1.5mm, Width: 5mm

2.0S Mini Reconstruction Locking Plate



REF	Holes	Length
021030004	4 holes	27mm
021030006	6 holes	40mm
021030008	8 holes	54mm
021030010	10 holes	67mm

Thickness: 1.5mm, Width: 5mm

2.0S Mini T Locking Plate



REF	Holes	Length
021040003	3 holes	26mm
021040005	5 holes	39mm
021040007	7 holes	53mm

Thickness: 1.5mm, Width: 5mm

2.0S Mini Y Locking Plate



REF	Holes	Length
021050003	3 holes	27mm
021050005	5 holes	41mm
021050007	7 holes	55mm

Thickness: 1.5mm, Width: 5mm

2.0S Mini T Reconstruction Locking Plate



REF	Holes	Length
021060003	3 holes	26mm
021060005	5 holes	39mm
021060007	7 holes	53mm

Thickness: 1.5mm, Width: 5mm

2.0S Mini Condylar Locking Plate



REF	Holes	Length
021070003	3 holes	27mm
021070005	5 holes	41mm
021070007	7 holes	55mm

Thickness: 1.5mm, Width: 5mm

Humeral Shaft Locking Plate



REF	Specification	Thickness	Width	Length
5100-0101	6 holes	3.6	13	92
5100-0102	7 holes	3.6	13	105
5100-0103	8 holes	3.6	13	118
5100-0104	9 holes	3.6	13	131
5100-0105	10 holes	3.6	13	144
5100-0106	12 holes	3.6	13	170
5100-0107	14 holes	3.6	13	196

Use 3.5 Locking Screw/ 3.5 Cortical Screw

1/3 Tubular Locking Plate



REF	Specification	Thickness	Width	Length
5100-0201	5 holes	2	10	71
5100-0202	6 holes	2	10	84
5100-0203	7 holes	2	10	97
5100-0204	8 holes	2	10	110
5100-0205	9 holes	2	10	123
5100-0206	10 holes	2	10	136
5100-0207	12 holes	2	10	162

Use 3.5 Locking Screw/ 3.5 Cortical Screw

Ulna & Radius Locking Plate



REF	Specification	Thickness	Width	Length
5100-0301	6 holes	3.2	11	92
5100-0302	7 holes	3.2	11	105
5100-0303	8 holes	3.2	11	118
5100-0304	9 holes	3.2	11	131
5100-0305	10 holes	3.2	11	144
5100-0306	12 holes	3.2	11	170

Use 3.5 Locking Screw/ 3.5 Cortical Screw

S-clavicle Locking Plate



REF	Specification	Thickness	Width	Length
5100-0401	5 holes L	3	10.5	54
5100-0402	6 holes L	3	10.5	67
5100-0403	7 holes L	3	10.5	81
5100-0404	8 holes L	3	10.5	94
5100-0405	9 holes L	3	10.5	107
5100-0406	10 holes L	3	10.5	131
5100-0407	5 holes R	3	10.5	54
5100-0408	6 holes R	3	10.5	67
5100-0409	7 holes R	3	10.5	81
5100-0410	8 holes R	3	10.5	94
5100-0411	9 holes R	3	10.5	107
5100-0412	10 holes R	3	10.5	131

Use 3.5 Locking Screw/ 3.5 Cortical Screw/ 4.0 Cancellous Screw

Distal Clavicular Locking Plate I



REF	Specification	Thickness	Width	Length
5100-0501	3 holes L	3	10	55
5100-0502	4 holes L	3	10	68
5100-0503	5 holes L	3	10	81
5100-0504	6 holes L	3	10	94
5100-0505	7 holes L	3	10	107
5100-0506	3 holes R	3	10	55
5100-0507	4 holes R	3	10	68
5100-0508	5 holes R	3	10	81
5100-0509	6 holes R	3	10	94
5100-0510	7 holes R	3	10	107

Use 2.7/3.5 Locking Screw/ 3.5 Cortical Screw/ 4.0 Cancellous Screw

Clavicle Hook Locking Plate



REF	Specification	Thickness	Width	Length
5100-0601	2 holes L	3	11.5	89
5100-0602	3 holes L	3	11.5	103
5100-0603	4 holes L	3	11.5	119
5100-0604	5 holes L	3	11.5	134
5100-0605	6 holes L	3	11.5	148
5100-0606	2 holes R	3	11.5	89
5100-0607	3 holes R	3	11.5	103
5100-0608	4 holes R	3	11.5	119
5100-0609	5 holes R	3	11.5	134
5100-0610	6 holes R	3	11.5	148

Use 3.5 Locking Screw/ 3.5 Cortical Screw/ 4.0 Cancellous Screw

Olecranon Locking Plate



REF	Specification	Thickness	Width	Length
5100-0701	3 holes L	2.5	11	107
5100-0702	4 holes L	2.5	11	120
5100-0703	6 holes L	2.5	11	146
5100-0704	8 holes L	2.5	11	172
5100-0705	10 holes L	2.5	11	198
5100-0706	3 holes R	2.5	11	107
5100-0707	4 holes R	2.5	11	120
5100-0708	6 holes R	2.5	11	146
5100-0709	8 holes R	2.5	11	172
5100-0710	10 holes R	2.5	11	198

Use 3.5 Locking Screw/ 3.5 Cortical Screw/ 4.0 Cancellous Screw

Distal Ulnar Locking Plate



REF	Specification	Thickness	Width	Length
5100-0801	6 holes	1.3	5.1	39
5100-0802	8 holes	1.3	5.1	49
5100-0803	10 holes	1.3	5.1	59

Use 2.4 Locking Screw/ 2.4 Cortical Screw

VA Distal Lateral Radius Locking Plate



REF	Specification	Thickness	Width	Length
5100-0901	5 holes	2	6.7	47
5100-0902	6 holes	2	6.7	55

Use 2.7 Locking Screw/ 2.7 Cortical Screw

VA Distal Medial Radius Locking Plate



REF	Specification	Thickness	Width	Length
5100-1001	4 holes L	2	7.2	41
5100-1002	5 holes L	2	7.2	48
5100-1003	6 holes L	2	7.2	55
5100-1004	4 holes R	2	7.2	41
5100-1005	5 holes R	2	7.2	48
5100-1006	6 holes R	2	7.2	55

Use 2.7 Locking Screw/ 2.7 Cortical screw

Distal Dorsal Radial Delta Locking Plate



REF	Specification	Thickness	Width	Length
S100-1101	3 holes L	1.8	8.6	65
S100-1102	4 holes L	1.8	8.6	73
S100-1103	5 holes L	1.8	8.6	82
S100-1104	7 holes L	1.8	8.6	99
S100-1105	9 holes L	1.8	8.6	116
S100-1106	3 holes R	1.8	8.6	65
S100-1107	4 holes R	1.8	8.6	73
S100-1108	5 holes R	1.8	8.6	82
S100-1109	7 holes R	1.8	8.6	99
S100-1110	9 holes R	1.8	8.6	116

Use 2.7 Locking Screw/ 2.7 Cortical screw

VA Distal Radius Locking Plate



REF	Specification	Thickness	Width	Length
S100-1201	3 holes L	2.2	8	60
S100-1202	4 holes L	2.2	8	69
S100-1203	5 holes L	2.2	8	77
S100-1204	6 holes L	2.2	8	86
S100-1205	7 holes L	2.2	8	94
S100-1206	3 holes R	2.2	8	60
S100-1207	4 holes R	2.2	8	69
S100-1208	5 holes R	2.2	8	77
S100-1209	6 holes R	2.2	8	86
S100-1210	7 holes R	2.2	8	94

Use 2.7 Locking Screw/ 2.7 Cortical Screw

Distal Volar Radial Locking Plate
(With Drill Guide)

REF	Specification	Thickness	Width	Length
S100-1301	3 holes L	2.5	9	49
S100-1302	4 holes L	2.5	9	58
S100-1303	5 holes L	2.5	9	66
S100-1304	7 holes L	2.5	9	83
S100-1305	9 holes L	2.5	9	99
S100-1306	3 holes R	2.5	9	49
S100-1307	4 holes R	2.5	9	58
S100-1308	5 holes R	2.5	9	66
S100-1309	7 holes R	2.5	9	83
S100-1310	9 holes R	2.5	9	99

Use 2.7 Locking Screw/ 2.7 Cortical Screw

Proximal Radius Locking Plate



REF	Specification	Thickness	Width	Length
S100-1401	3 holes L	1.8	8.7	53
S100-1402	4 holes L	1.8	8.7	63
S100-1403	5 holes L	1.8	8.7	72
S100-1404	3 holes R	1.8	8.7	53
S100-1405	4 holes R	1.8	8.7	63
S100-1406	5 holes R	1.8	8.7	72

Use 2.4 Locking Screw/ 2.4 Cortical Screw

Proximal Humeral Locking Plate



REF	Specification	Thickness	Width	Length
S100-1501	3 holes	4	12	90
S100-1502	4 holes	4	12	102
S100-1503	5 holes	4	12	114
S100-1504	6 holes	4	12	126
S100-1505	7 holes	4	12	138
S100-1506	8 holes	4	12	150
S100-1507	10 holes	4	12	174
S100-1508	12 holes	4	12	198

Use 3.5 Locking Screw/ 3.5 Cortical Screw/ 4.0 Cancellous Screw

Proximal Humeral Greater Tuberosity
Locking Plate

REF	Specification	Thickness	Width	Length
S100-1601	5 holes L	1.5	13	44
S100-1602	5 holes R	1.5	13	44

Use 2.7/3.5 Locking Screw, 2.7/3.5 Cortical Screw/ 4.0 Cancellous Screw

Distal Lateral Humeral Locking Plate



REF	Specification	Thickness	Width	Length
S100-1701	4 holes L	3.2	12	86
S100-1702	6 holes L	3.2	12	112
S100-1703	8 holes L	3.2	12	138
S100-1704	10 holes L	3.2	12	164
S100-1705	12 holes L	3.2	12	190
S100-1706	4 holes R	3.2	12	86
S100-1707	6 holes R	3.2	12	112
S100-1708	8 holes R	3.2	12	138
S100-1709	10 holes R	3.2	12	164
S100-1710	12 holes R	3.2	12	190

Use 2.7/3.5 Locking Screw/ 3.5 Cortical Screw

Distal Medial Humeral Locking Plate



REF	Specification	Thickness	Width	Length
S100-1801	4 holes L	3	11.5	69
S100-1802	6 holes L	3	11.5	95
S100-1803	8 holes L	3	11.5	121
S100-1804	10 holes L	3	11.5	147
S100-1805	12 holes L	3	11.5	173
S100-1806	4 holes R	3	11.5	69
S100-1807	6 holes R	3	11.5	95
S100-1808	8 holes R	3	11.5	121
S100-1809	10 holes R	3	11.5	147
S100-1810	12 holes R	3	11.5	173

Use 2.7/3.5 Locking Screw/ 3.5 Cortical Screw

Distal Fibular Locking Plate



REF	Specification	Thickness	Width	Length
S100-1901	3 holes L	2.5	10	81
S100-1902	4 holes L	2.5	10	93
S100-1903	5 holes L	2.5	10	105
S100-1904	6 holes L	2.5	10	117
S100-1905	7 holes L	2.5	10	129
S100-1906	8 holes L	2.5	10	141
S100-1907	3 holes R	2.5	10	81
S100-1908	4 holes R	2.5	10	93
S100-1909	5 holes R	2.5	10	105
S100-1910	6 holes R	2.5	10	117
S100-1911	7 holes R	2.5	10	129
S100-1912	8 holes R	2.5	10	141

Use 3.5 Locking Screw/ 3.5 Cortical Screw

Distal Clavicular Locking Plate



REF	Holes	Length
01.1158-01	3 holes L	61mm
01.1158-02	4 holes L	71mm
01.1158-03	5 holes L	81mm
01.1158-11	3 holes R	61mm
01.1158-12	4 holes R	71mm
01.1158-13	5 holes R	81mm

Use 2.7/3.5 Locking Screw/3.5 Cortical Screw

Olecranon Locking Plate II



REF	Holes	Length
01.1062-01	3 holes	49mm
01.1062-02	4 holes	63mm
01.1062-03	5 holes	77mm
01.1062-04	6 holes	91mm
01.1062-05	7 holes	105mm
01.1062-06	8 holes	119mm

Use 3.5 Locking Screw/3.5 Cortical Screw

Distal Dorsal Radial Delta Locking Plate



REF	Holes	Length
01.1065-01	3 holes L	46mm
01.1065-02	4 holes L	56mm
01.1065-03	5 holes L	66mm
01.1065-04	6 holes L	76mm
01.1065-11	3 holes R	46mm
01.1065-12	4 holes R	56mm
01.1065-13	5 holes R	66mm
01.1065-14	6 holes R	76mm

Use 2.7 Locking Screw/2.7 Cortical Screw

Distal Humeral Extra-articular Locking Plate



REF	Holes	Length
01.1107-01	4 holes L	114mm
01.1107-02	5 holes L	128mm
01.1107-03	6 holes L	142mm
01.1107-04	7 holes L	156mm
01.1107-05	8 holes L	170mm
01.1107-06	9 holes L	184mm
01.1107-11	4 holes R	114mm
01.1107-12	5 holes R	128mm
01.1107-13	6 holes R	142mm
01.1107-14	7 holes R	156mm
01.1107-15	8 holes R	170mm
01.1107-16	9 holes R	184mm

Use 3.5 Locking Screw/3.5 Cortical Screw

Distal Fibular Locking Plate-II



REF	Holes	Length
01.1222-01	3 holes L	84mm
01.1222-02	4 holes L	97mm
01.1222-03	5 holes L	110mm
01.1222-04	6 holes L	123mm
01.1222-05	7 holes L	136mm
01.1222-06	9 holes L	162mm
01.1222-07	11 holes L	188mm
01.1222-11	3 holes R	84mm
01.1222-12	4 holes R	97mm
01.1222-13	5 holes R	110mm
01.1222-14	6 holes R	123mm
01.1222-15	7 holes R	136mm
01.1222-16	9 holes R	162mm
01.1222-17	11 holes R	188mm

Use 2.7/3.5 Locking Screw/2.7/3.5 Cortical Screw

Proximal Lateral Tibial Osteotomy Locking Plate



REF	Holes	Length
01.1223-01	3 holes L	104mm
01.1223-11	3 holes R	104mm

Use 5.0 Locking Screw/4.5 Cortical Screw

Distal Medial Femoral Osteotomy Locking Plate



REF	Holes	Length
01.1275-01	3 holes L	121mm
01.1275-02	4 holes L	101mm
01.1275-11	3 holes R	121mm
01.1275-12	4 holes R	101mm

Use 5.0 Locking Screw/4.5 Cortical Screw

Distal Medial Femoral Osteotomy Locking Plate



REF	Holes	Length
5100-5403	5 holes L	161mm
5100-5404	5 holes R	161mm

Use 5.0 Locking Screw/4.5 Cortical Screw

Lateral Ankle Fusion Locking Plate



REF	Holes	Length
01.1224-01	4 holes L	98mm
01.1224-11	4 holes R	98mm

Use 3.5 Locking Screw/3.5 Cortical Screw

Tibia Locking Plate Narrow



REF	Specification	Thickness	Width	Length
5100-2201	7 holes	5.0	15	139
5100-2202	8 holes	5.0	15	157
5100-2203	9 holes	5.0	15	175
5100-2204	10 holes	5.0	15	193
5100-2205	12 holes	5.0	15	229
5100-2206	14 holes	5.0	15	265

Use 5.0 Locking Screw/ 4.5 Cortical Screw

Proximal Medial Tibial Osteotomy Locking Plate



REF	Specification	Thickness	Width	Length
5100-2301	5 holes	2.8	16	115

Use 5.0 Locking Screw/ 4.5 Cortical Screw

Lateral Tibial Head Buttress Locking Plate



REF	Specification	Thickness	Width	Length
5100-2401	5 holes L	4.6	15	144
5100-2402	7 holes L	4.6	15	182
5100-2403	9 holes L	4.6	15	220
5100-2404	11 holes L	4.6	15	258
5100-2405	13 holes L	4.6	15	296
5100-2406	5 holes R	4.6	15	144
5100-2407	7 holes R	4.6	15	182
5100-2408	9 holes R	4.6	15	220
5100-2409	11 holes R	4.6	15	258
5100-2410	13 holes R	4.6	15	296

Use 5.0 Locking Screw/ 4.5 Cortical Screw

Proximal Lateral Tibial Locking Plate-I



REF	Specification	Thickness	Width	Length
5100-2501	3 holes L	4.6	14	117
5100-2502	5 holes L	4.6	14	155
5100-2503	7 holes L	4.6	14	193
5100-2504	9 holes L	4.6	14	231
5100-2505	11 holes L	4.6	14	269
5100-2506	3 holes R	4.6	14	117
5100-2507	5 holes R	4.6	14	155
5100-2508	7 holes R	4.6	14	193
5100-2509	9 holes R	4.6	14	231
5100-2510	11 holes R	4.6	14	269

Use 3.5/5.0 Locking Screw/ 4.5 Cortical Screw

Proximal Lateral Tibial Locking Plate-II



REF	Specification	Thickness	Width	Length
5100-2601	5 holes L	4.2	12.5	97
5100-2602	7 holes L	4.2	12.5	124
5100-2603	9 holes L	4.2	12.5	151
5100-2604	11 holes L	4.2	12.5	178
5100-2605	13 holes L	4.2	12.5	205
5100-2606	5 holes R	4.2	12.5	97
5100-2607	7 holes R	4.2	12.5	124
5100-2608	9 holes R	4.2	12.5	151
5100-2609	11 holes R	4.2	12.5	178
5100-2610	13 holes R	4.2	12.5	205

Use 3.5 Locking Screw/ 3.5 Cortical Screw

Proximal Medial Tibial Locking Plate



REF	Specification	Thickness	Width	Length
5100-2701	4 holes L	4.2	13	83
5100-2702	6 holes L	4.2	13	109
5100-2703	8 holes L	4.2	13	135
5100-2704	10 holes L	4.2	13	161
5100-2705	12 holes L	4.2	13	187
5100-2706	4 holes R	4.2	13	83
5100-2707	6 holes R	4.2	13	109
5100-2708	8 holes R	4.2	13	135
5100-2709	10 holes R	4.2	13	161
5100-2710	12 holes R	4.2	13	187

Use 3.5 Locking Screw/ 3.5 Cortical Screw

Distal Lateral Tibial Locking Plate-I



REF	Specification	Thickness	Width	Length
5100-2801	5 holes L	3.6	16.5	122
5100-2802	7 holes L	3.6	16.5	154
5100-2803	9 holes L	3.6	16.5	186
5100-2804	11 holes L	3.6	16.5	218
5100-2805	13 holes L	3.6	16.5	250
5100-2806	5 holes R	3.6	16.5	122
5100-2807	7 holes R	3.6	16.5	154
5100-2808	9 holes R	3.6	16.5	186
5100-2809	11 holes R	3.6	16.5	218
5100-2810	13 holes R	3.6	16.5	250

Use 5.0 Locking Screw/ 4.5 Cortical Screw

Distal Lateral Tibial Locking Plate-II



REF	Specification	Thickness	Width	Length
5100-2901	5 holes L	4.2	13	136
5100-2902	7 holes L	4.2	13	168
5100-2903	9 holes L	4.2	13	200
5100-2904	11 holes L	4.2	13	232
5100-2905	5 holes R	4.2	13	136
5100-2906	7 holes R	4.2	13	168
5100-2907	9 holes R	4.2	13	200
5100-2908	11 holes R	4.2	13	232

Use 3.5 Locking Screw/ 3.5 Cortical Screw

Distal Medial Tibial Locking Plate-I



REF	Specification	Thickness	Width	Length
5100-3001	5 holes L	4.2	14	147
5100-3002	7 holes L	4.2	14	179
5100-3003	9 holes L	4.2	14	211
5100-3004	11 holes L	4.2	14	243
5100-3005	13 holes L	4.2	14	275
5100-3006	5 holes R	4.2	14	147
5100-3007	7 holes R	4.2	14	179
5100-3008	9 holes R	4.2	14	211
5100-3009	11 holes R	4.2	14	243
5100-3010	13 holes R	4.2	14	275

Use 3.5 Locking Screw/ 3.5 Cortical Screw

Distal Medial Tibial Locking Plate-II



REF	Specification	Thickness	Width	Length
S100-3101	5 holes L	4.2	15.5	125
S100-3102	7 holes L	4.2	15.5	151
S100-3103	9 holes L	4.2	15.5	177
S100-3104	11 holes L	4.2	15.5	203
S100-3105	5 holes R	4.2	15.5	125
S100-3106	7 holes R	4.2	15.5	151
S100-3107	9 holes R	4.2	15.5	177
S100-3108	11 holes R	4.2	15.5	203

Use 3.5 Locking Screw/ 3.5 Cortical Screw

Curved Femur Locking Compression Plate
(Broad)

REF	Specification	Thickness	Width	Length
S100-3201	8 holes	6.0	18	158
S100-3202	9 holes	6.0	18	176
S100-3203	10 holes	6.0	18	194
S100-3204	12 holes	6.0	18	230
S100-3205	14 holes	6.0	18	266
S100-3206	16 holes	6.0	18	302

Use 5.0 Locking Screw/ 4.5 Cortical Screw

Straight Reconstruction Locking Plate



REF	Specification	Thickness	Width	Length
S100-3301	5 holes	3.2	11	66
S100-3302	6 holes	3.2	11	79
S100-3303	7 holes	3.2	11	92
S100-3304	8 holes	3.2	11	105
S100-3305	9 holes	3.2	11	118
S100-3306	10 holes	3.2	11	131
S100-3307	12 holes	3.2	11	157

Use 3.5 Locking Screw/ 3.5 Cortical Screw

Patella Mesh Locking Plate



REF	Specification	Thickness	Width	Length
S100-3401	16 holes	Small	30	38
S100-3402	16 holes	Medium	33	42
S100-3403	16 holes	Large	36	46

Use 2.7 Locking Screw

Distal Femoral Locking Plate



REF	Specification	Thickness	Width	Length
S100-3501	5 holes L	6.0	17.5	164
S100-3502	7 holes L	6.0	17.5	204
S100-3503	9 holes L	6.0	17.5	244
S100-3504	11 holes L	6.0	17.5	284
S100-3505	13 holes L	6.0	17.5	324
S100-3506	5 holes R	6.0	17.5	164
S100-3507	7 holes R	6.0	17.5	204
S100-3508	9 holes R	6.0	17.5	244
S100-3509	11 holes R	6.0	17.5	284
S100-3510	13 holes R	6.0	17.5	324

Use 5.0 Locking Screw/ 4.5 Cortical Screw

Proximal Femoral Locking Plate-I



REF	Specification	Thickness	Width	Length
S100-3601	5 holes L	6.0	18	153
S100-3602	7 holes L	6.0	18	189
S100-3603	9 holes L	6.0	18	225
S100-3604	11 holes L	6.0	18	261
S100-3605	13 holes L	6.0	18	297
S100-3606	5 holes R	6.0	18	153
S100-3607	7 holes R	6.0	18	189
S100-3608	9 holes R	6.0	18	225
S100-3609	11 holes R	6.0	18	261
S100-3610	13 holes R	6.0	18	297

Use 5.0 Locking Screw/ 4.5 Cortical Screw

Proximal Femur Locking Plate-II



REF	Specification	Thickness	Width	Length
S100-3701	3 holes L	6.5	19	133
S100-3702	4 holes L	6.5	19	151
S100-3703	5 holes L	6.5	19	169
S100-3704	7 holes L	6.5	19	205
S100-3705	9 holes L	6.5	19	241
S100-3706	11 holes L	6.5	19	277
S100-3707	3 holes R	6.5	19	133
S100-3708	4 holes R	6.5	19	151
S100-3709	5 holes R	6.5	19	169
S100-3710	7 holes R	6.5	19	205
S100-3711	9 holes R	6.5	19	241
S100-3712	11 holes R	6.5	19	277

Use 5.0 Locking Screw/ 4.5 Cortical Screw/ 6.5 Cannulated Full-threaded Locking Screw

Calcaneus Locking Plate



REF	Specification	Thickness	Width	Length
S100-3801	Small Right	2	34	60
S100-3802	Small Left	2	34	60
S100-3803	Medium Right	2	34.5	67
S100-3804	Medium Left	2	34.5	67
S100-3805	Large Right	2	35	73
S100-3806	Large Left	2	35	73

Use 3.5 Locking Screw/ 4.0 Cortical Screw

Acromion of Scapula Locking Plate



REF	Holes	Length
01.1154-01	9 holes L	79mm
01.1154-11	9 holes R	79mm

Use 2.7 Locking Screw

Scapula Medial Locking Plate



REF	Holes	Length
01.1157-01	10 holes L	72mm
01.1157-11	10 holes R	72mm

Use 2.7 Locking Screw

Scapula Lateral Locking Plate



REF	Holes	Length
01.1155-01	6 holes L	89mm
01.1155-11	6 holes R	89mm

Use 2.7 Locking Screw

Glenoid Cavity Locking Plate



REF	Holes	Thickness
01.1156-01	4 holes L	1.5mm
01.1156-11	4 holes R	1.5mm

Use 2.7 Locking Screw

Wrist Fusion Locking Plate



REF	Holes	Length
01.1067-01	9 holes	117mm
01.1067-11	8 holes	112mm

Use 2.7/3.5 Locking Screw/2.7/3.5 Cortical Screw

Distal Tibial T-Locking Plate



REF	Holes	Length
5100-6303	3 holes	56.5mm
5100-6304	4 holes	73.5mm
5100-6305	6 holes	93.5mm
5100-6306	8 holes	113.5mm

Use 3.5 Locking Screw/3.5 Cortical Screw

MD Distal Tibial Locking Plate



REF	Holes	Length
5100-6307	3 holes	53mm
5100-6308	4 holes	65mm

Use 3.5 Locking Screw/3.5 Cortical Screw

Distal Tibial Medial Locking Plate



REF	Holes	Length
5100-6301	4 holes	88mm
5100-6302	6 holes	112mm

Use 3.5 Locking Screw/3.5 Cortical Screw

MD Distal Tibial Locking Plate



REF	Holes	Length
01.1220-01	3 holes L	70mm
01.1220-02	4 holes L	70mm
01.1220-11	3 holes R	93mm
01.1220-12	4 holes R	93mm

Use 3.5 Locking Screw/3.5 Cortical Screw

Proximal Femoral Hook Locking Plate with Cable



REF	Holes	Length
5100-4701	4 holes L	233mm
5100-4702	6 holes L	258mm
5100-4703	8 holes L	283mm
5100-4704	10 holes L	308mm
5100-4705	12 holes L	333mm
5100-4706	4 holes R	233mm
5100-4707	6 holes R	258mm
5100-4708	8 holes R	283mm
5100-4709	10 holes R	308mm
5100-4710	12 holes R	333mm

Use 5.0 Locking Screw/ 4.5 Cortical Screw

Proximal Femoral Hook Locking Plate



REF	Holes	Length
5100-5301	3 holes L	157.5mm
5100-5302	4 holes L	175.5mm
5100-5303	5 holes L	193.5mm
5100-5304	7 holes L	229.5mm
5100-5305	9 holes L	265.5mm
5100-5306	11 holes L	301.5mm
5100-5307	13 holes L	337.5mm
5100-5308	3 holes R	157.5mm
5100-5309	4 holes R	175.5mm
5100-5310	5 holes R	193.5mm
5100-5311	7 holes R	229.5mm
5100-5312	9 holes R	265.5mm
5100-5313	11 holes R	301.5mm
5100-5314	13 holes R	337.5mm

Use 5.0 Locking Screw/ 4.5 Cortical Screw/ 6.5 Cannula Locking Screw

Cable Positioning Pin



REF	Specification
5100-4803	3.5mm
5100-4804	5.0mm

Subtalar Screw (Tarsal Tunnel Approach)



REF	Specification	REF	Specification
3105355021	Φ5.5×21mm	3105385021	Φ8.5×21
3105365021	Φ6.5×21mm	3105395021	Φ9.5×21
3105373021	Φ7.3×21mm		

Subtalar Screw (Tarsal Sinus Approach)



REF	Specification	REF	Specification
3152385015	Φ8.5*15mm	3152311020	Φ11*20mm
3152310015	Φ10*15mm	3152312015	Φ12*15mm
3152311015	Φ11*15mm	3152312020	Φ12*20mm

Distal Tibial Medial Locking Plate



REF	Specification	Length
01.4070-01	Small Right	56.3mm
01.4070-02	Small Left	58.9mm
01.4070-03	Medium Right	67.8mm
01.4070-11	Medium Left	56.3mm
01.4070-12	Large Right	58.9mm
01.4070-13	Large Left	67.8mm

Use 3.5 Locking Screw

MD Posterior Tuberosity Calcaneal Locking Plate



REF	Specification	Length
01.4071-01	Small Right	39.8mm
01.4071-02	Small Left	49.3mm
01.4071-11	Large Right	39.8mm
01.4071-12	Large Left	49.3mm

Use 3.5 Locking Screw

MD Anterior Process Calcaneal Locking Plate



REF	Specification	Length
01.4072-01	Small Right	46mm
01.4072-02	Small Left	50.8mm
01.4072-11	Large Right	46mm
01.4072-12	Large Left	50.8mm

Use 3.5 Locking Screw

MT Fifth Metatarsal Locking Plate



REF	Specification	Length
01.4073-01	6 holes L	48mm
01.4073-11	6 holes R	48mm

Use 2.7 Locking Screw/ 2.7 Cortical Screw

MT First Metatarsocuneiform Fusion Locking Plate



REF	Specification	Length
01.4074-01	8 holes	60mm

Use 2.7 Locking Screw/ 2.7 Cortical Screw

MT L-shape Metatarsal Locking Plate



REF	Specification	Length
01.4075-01	4 holes L	45.1mm
01.4075-02	5 holes L	51.6mm
01.4075-11	4 holes R	45.1mm
01.4075-12	5 holes R	51.6mm

Use 2.7 Locking Screw/ 2.7 Cortical Screw

MT L-shape Metatarsal Locking Plate



MT Straight Fusion Locking Plate



MT First Metatarsophalangeal Fusion Locking Plate



MT X-Locking Plate



MT Cuboid Locking Plate



MT Astragalus Locking Plate



MT Navicular Locking Plate



REF	Specification	Length
01.4076-01	2 holes	27.4mm
01.4076-02	3 holes	33.9mm
01.4076-03	4 holes	40.4mm
01.4076-04	5 holes	47mm
01.4076-05	6 holes	53.7mm

Use 2.7 Locking Screw/ 2.7 Cortical Screw

REF	Specification	Length
01.4077-01	2 holes	18.5mm
01.4077-02	4 holes	35.6mm
01.4077-03	3 holes	30mm

Use 2.7 Locking Screw/ 2.7 Cortical Screw

REF	Specification	Length
01.4078-01	Small Right	55.3mm
01.4078-02	Small Left	55.3mm
01.4078-11	Large Right	60mm
01.4078-12	Large Left	60mm

Use 2.7 Locking Screw/ 2.7 Cortical Screw

REF	Specification	Length
01.4079-01	24mm	24mm
01.4079-02	28.5mm	28.5mm

Use 2.7 Locking Screw

REF	Specification	Length
01.4080-01	9 holes L	29.5mm
01.4080-11	9 holes R	29.5mm

Use 2.7 Locking Screw

REF	Specification	Length
01.4081-01	5 holes	43mm

Use 2.7 Locking Screw

REF	Specification	Length
01.4082-01	9 holes L	46mm
01.4082-11	9 holes R	46mm

Use 2.7 Locking Screw

1.5mm Cortex Screw
(T4 Stardrive) (Self-tapping)



REF	Length
030310006	1.5*6mm
030310008	1.5*8mm
030310010	1.5*10mm
030310012	1.5*12mm
030310014	1.5*14mm
030310016	1.5*16mm
030310018	1.5*18mm
030310020	1.5*20mm

1.5mm Locking Screw
(T4 Stardrive) (Self-tapping)



REF	Length
030300006	1.5*6mm
030300008	1.5*8mm
030300010	1.5*10mm
030300012	1.5*12mm
030300014	1.5*14mm
030300016	1.5*16mm
030300018	1.5*18mm
030300020	1.5*20mm

2.0mm Cortex Screw
(T6 Stardrive) (Self-tapping)



REF	Length
030330006	2.0*6mm
030330008	2.0*8mm
030330010	2.0*10mm
030330012	2.0*12mm
030330014	2.0*14mm
030330016	2.0*16mm
030330018	2.0*18mm
030330020	2.0*20mm
030330022	2.0*22mm
030330024	2.0*24mm
030330026	2.0*26mm
030330028	2.0*28mm
030330030	2.0*30mm

2.0mm Locking Screw
(T6 Stardrive) (Self-tapping)



REF	Length
030320006	2.0*6mm
030320008	2.0*8mm
030320010	2.0*10mm
030320012	2.0*12mm
030320014	2.0*14mm
030320016	2.0*16mm
030320018	2.0*18mm
030320020	2.0*20mm
030320022	2.0*22mm
030320024	2.0*24mm
030320026	2.0*26mm
030320028	2.0*28mm
030320030	2.0*30mm

2.4mm Cortex Screw
(T8 Stardrive) (Self-tapping)



REF	Length
030350006	2.4*6mm
030350008	2.4*8mm
030350010	2.4*10mm
030350012	2.4*12mm
030350014	2.4*14mm
030350016	2.4*16mm
030350018	2.4*18mm
030350020	2.4*20mm
030350022	2.4*22mm
030350024	2.4*24mm
030350026	2.4*26mm
030350028	2.4*28mm
030350030	2.4*30mm

2.4mm Locking Screw
(T8 Stardrive) (Self-tapping)



REF	Length
040240006	2.4*6mm
040240008	2.4*8mm
040240010	2.4*10mm
040240012	2.4*12mm
040240014	2.4*14mm
040240016	2.4*16mm
040240018	2.4*18mm
040240020	2.4*20mm
040240022	2.4*22mm
040240024	2.4*24mm
040240026	2.4*26mm
040240028	2.4*28mm
040240030	2.4*30mm

2.7mm Cortex Screw
(T8 Stardrive) (Self-tapping)



REF	Length
030390010	2.7*10mm
030390012	2.7*12mm
030390014	2.7*14mm
030390016	2.7*16mm
030390018	2.7*18mm
030390020	2.7*20mm
030390022	2.7*22mm
030390024	2.7*24mm
030390026	2.7*26mm
030390028	2.7*28mm
030390030	2.7*30mm

2.7mm Locking Screw
(T8 Stardrive) (Self-tapping)



REF	Length	REF	Length
030340010	2.7*10mm	030340026	2.7*26mm
030340012	2.7*12mm	030340028	2.7*28mm
030340014	2.7*14mm	030340030	2.7*30mm
030340016	2.7*16mm	030340032	2.7*32mm
030340018	2.7*18mm	030340034	2.7*34mm
030340020	2.7*20mm	030340036	2.7*36mm
030340022	2.7*22mm	030340038	2.7*38mm
030340024	2.7*24mm	030340040	2.7*40mm

3.5mm Locking Screw (Stardrive)



REF	Specification
5100-3901	3.5*12mm
5100-3902	3.5*14mm
5100-3903	3.5*16mm
5100-3904	3.5*18mm
5100-3905	3.5*20mm
5100-3906	3.5*22mm
5100-3907	3.5*24mm
5100-3908	3.5*26mm
5100-3909	3.5*28mm
5100-3910	3.5*30mm
5100-3911	3.5*32mm
5100-3912	3.5*34mm
5100-3913	3.5*36mm
5100-3914	3.5*38mm
5100-3915	3.5*40mm
5100-3916	3.5*42mm

REF	Specification
5100-3917	3.5*44mm
5100-3918	3.5*46mm
5100-3919	3.5*48mm
5100-3920	3.5*50mm
5100-3921	3.5*52mm
5100-3922	3.5*54mm
5100-3923	3.5*56mm
5100-3924	3.5*58mm
5100-3925	3.5*60mm
5100-3926	3.5*65mm
5100-3927	3.5*70mm
5100-3928	3.5*75mm
5100-3929	3.5*80mm
5100-3930	3.5*85mm
5100-3931	3.5*90mm
5100-3932	3.5*95mm

5.0mm Locking Screw (Stardrive)



REF	Specification
5100-4001	5.0*22mm
5100-4002	5.0*24mm
5100-4003	5.0*26mm
5100-4004	5.0*28mm
5100-4005	5.0*30mm
5100-4006	5.0*32mm
5100-4007	5.0*34mm
5100-4008	5.0*36mm
5100-4009	5.0*38mm
5100-4010	5.0*40mm
5100-4011	5.0*42mm
5100-4012	5.0*44mm
5100-4013	5.0*46mm
5100-4014	5.0*48mm

REF	Specification
5100-4015	5.0*50mm
5100-4016	5.0*52mm
5100-4017	5.0*54mm
5100-4018	5.0*56mm
5100-4019	5.0*58mm
5100-4020	5.0*60mm
5100-4021	5.0*65mm
5100-4022	5.0*70mm
5100-4023	5.0*75mm
5100-4024	5.0*80mm
5100-4025	5.0*85mm
5100-4026	5.0*90mm
5100-4027	5.0*95mm

3.5mm Cortical Screw (Stardrive)



REF	Specification
5100-4101	3.5*12mm
5100-4102	3.5*14mm
5100-4103	3.5*16mm
5100-4104	3.5*18mm
5100-4105	3.5*20mm
5100-4106	3.5*22mm
5100-4107	3.5*24mm
5100-4108	3.5*26mm
5100-4109	3.5*28mm
5100-4110	3.5*30mm
5100-4111	3.5*32mm

REF	Specification
5100-4112	3.5*34mm
5100-4113	3.5*36mm
5100-4114	3.5*38mm
5100-4115	3.5*40mm
5100-4116	3.5*42mm
5100-4117	3.5*44mm
5100-4118	3.5*46mm
5100-4119	3.5*48mm
5100-4120	3.5*50mm
5100-4121	3.5*55mm
5100-4122	3.5*60mm

4.5mm Cortical Screw (Stardrive)



REF	Specification
5100-4201	4.5*22mm
5100-4202	4.5*24mm
5100-4203	4.5*26mm
5100-4204	4.5*28mm
5100-4205	4.5*30mm
5100-4206	4.5*32mm
5100-4207	4.5*34mm
5100-4208	4.5*36mm
5100-4209	4.5*38mm
5100-4210	4.5*40mm

REF	Specification
5100-4211	4.5*42mm
5100-4212	4.5*44mm
5100-4213	4.5*46mm
5100-4214	4.5*48mm
5100-4215	4.5*50mm
5100-4216	4.5*52mm
5100-4217	4.5*54mm
5100-4218	4.5*56mm
5100-4219	4.5*58mm
5100-4220	4.5*60mm

4.0mm Cancellous Screw (Stardrive)



REF	Specification
5100-4301	4.0*12mm
5100-4302	4.0*14mm
5100-4303	4.0*16mm
5100-4304	4.0*18mm
5100-4305	4.0*20mm
5100-4306	4.0*22mm
5100-4307	4.0*24mm
5100-4308	4.0*26mm
5100-4309	4.0*28mm
5100-4310	4.0*30mm
5100-4311	4.0*32mm
5100-4312	4.0*34mm
5100-4313	4.0*36mm
5100-4314	4.0*38mm
5100-4315	4.0*40mm
5100-4316	4.0*42mm
5100-4317	4.0*44mm
5100-4318	4.0*46mm
5100-4319	4.0*48mm
5100-4320	4.0*50mm

6.5 Cannulated Full-threaded Locking Screw (Stardrive)



REF	Specification
5100-4401	6.5*50mm
5100-4402	6.5*55mm
5100-4403	6.5*60mm
5100-4404	6.5*65mm
5100-4405	6.5*70mm
5100-4406	6.5*75mm
5100-4407	6.5*80mm
5100-4408	6.5*85mm
5100-4409	6.5*90mm
5100-4410	6.5*95mm
5100-4411	6.5*100mm
5100-4412	6.5*105mm
5100-4413	6.5*110mm

Rib Reconstruction Locking Plate



REF	Specification
5100-2001	15 holes L
5100-2002	18 holes L
5100-2003	15 holes R
5100-2004	18 holes R

Rib Plate Locking Screw (Self-tapping)



REF	Specification
5100-2101	2.9*8mm
5100-2102	2.9*10mm
5100-2103	2.9*12mm
5100-2104	2.9*14mm
5100-2105	2.9*16mm
5100-2106	2.9*18mm
5100-2107	2.9*20mm

FNS (Femoral Neck System)



Plate (130° CCD angle)



REF	Holes
3300-0101	1 hole
3300-0102	2 holes

Bolts+Antirotation Screws



REF	Length	REF	Length
3300-0201	75mm	3300-0206	100mm
3300-0202	80mm	3300-0207	105mm
3300-0203	85mm	3300-0208	110mm
3300-0204	90mm	3300-0209	115mm
3300-0205	95mm		

5.0mm Locking Screws (Self-tapping)



REF	Length	REF	Length
5100-4005	30mm	5100-4013	46mm
5100-4006	32mm	5100-4014	48mm
5100-4007	34mm	5100-4015	50mm
5100-4008	36mm	5100-4016	52mm
5100-4009	38mm	5100-4017	54mm
5100-4010	40mm	5100-4018	56mm
5100-4011	42mm	5100-4019	58mm
5100-4012	44mm	5100-4020	60mm

Sternum Locking Plate with Clip-I



REF	Thickness	Length
5100-4508	2.3mm	20mm

Sternum Locking Plate with Clip-II



REF	Thickness	Length
5100-4509	2.3mm	29mm

Sternum Locking Plate with Clip-III



REF	Thickness	Length
5100-4510	2.3mm	20mm

Sternum Locking Plate with Clip-IV



REF	Thickness	Length
5100-4511	2.3mm	24mm

Sternum Locking Plate with Clip-V



REF	Thickness	Length
5100-4512	2.3mm	33mm

Straight Sternum Locking Plate with Clip-I



REF	Thickness	Length
5100-4513	2.3mm	105mm

Straight Sternum Locking Plate with Clip-II




REF	Thickness	Length
5100-4514	2.3mm	169mm

Mini Fragment Locking Plate Instrument Set 3200-10






NO.	REF	Description	Qty.
1	3200-1001	1.5mm Plate Tray	1
2	3200-1002	2.0mm Plate Tray	1
3	3200-1003	2.4mm Plate Tray	1
4	3200-1004	1.5mm Plate Holder	2
5	3200-1005	1.1mm AO Drill Bit	1
6	3200-1006	T5 Stardriver	1
7	3200-1007	3.0mm Countersink	1
8	3200-1008	1.5mm AO Tap	1
9	3200-1009	1.1/1.5 Guider	1
10	3200-1010	1.1mm Locking Sleeve	1
11	3200-1011	1.1mm Locking Sleeve	1
12	3200-1012	2.0mm Plate Holder	2
13	3200-1013	1.5mm AO Drill Bit	1
14	3200-1014	T6 Stardriver	1
15	3200-1015	2.0mm AO Tap	1
16	3200-1016	4.0mm Countersink	1
17	3200-1017	1.5/2.0 Guider	1
18	3200-1018	1.5mm Locking Sleeve	1
19	3200-1019	1.5mm Locking Sleeve	1
20	3200-1020	2.4mm Plate Holder	2
21	3200-1021	2.0mm AO Drill Bit	1
22	3200-1022	T8 Stardriver	1
23	3200-1023	4.5mm Countersink	1
24	3200-1024	2.4mm AO Tap	1
25	3200-1025	2.0/2.4 Guider	1
26	3200-1026	2.0mm Locking Sleeve	1
27	3200-1027	2.0mm Locking Sleeve	1
28	3200-1028	Bending Forcep	1
29	3200-1029	Plate Cutter	1
30	3200-1030	0.4 N.m Torque Wrench	1
31	3200-1031	0.8 N.m Torque Wrench	1
32	3200-1032	Self-locking Reduction Forcep	1
33	3200-1033	Screwholding Forcep	1
34	3200-1034	AO Quick Coupling Handle	1
35	3200-1035	Tip Reduction Forcep	1
36	3200-1036	Depth Gauge 0-30mm	1
37	3200-1037	Depth Gauge 0-40mm	1
38	3200-1038	Periosteal Elevator Flat	1
39	3200-1039	Periosteal Elevator Round	1
40	3200-1040	Aluminium Box	1

2.7mm Mini Fragment Locking Plate Instrument Set 3200-11





NO.	REF	Description	Qty.
1	3200-1101	2.7mm Plate Tray	1
2	3200-1102	2.0/2.7 Guider	1
3	3200-1103	4.5mm Countersink	1
4	3200-1104	2.7mm AO Tap	1
5	3200-1105	T8 Stardriver	1
6	3200-1106	2.0mm AO Drill Bit	1
7	3200-1107	2.7mm Plate Holder	2
8	3200-1108	2.7mm Locking Sleeve	1
9	3200-1109	2.7mm Locking Sleeve	1

Small Fragment Locking Plate Instrument Set (AO) 3200-13






NO.	REF	Description	Qty.
1	3200-1301	Hex Key	1
2	3200-1302	Locking Sleeve 3.2mm for Thin Plates	1
3	3200-1303	Locking Sleeve 3.2mm for Thin Plates	1
4	3200-1304	Wire Sleeve 1.5mm	1
5	3200-1305	Locking Sleeve 3.2mm for Thick Plates	1
6	3200-1306	Locking Sleeve 3.2mm for Thick Plates	1
7	3200-1307	Locking Sleeve 3.2mm for Thick Plates	1
8	3200-1308	Transition Interface	1
9	3200-1309	Wire Sleeve 1.5mm	1
10	3200-1310	AO Drill Bit Ø2.5*200	1
11	3200-1311	AO Drill Bit Ø2.5*200	1
12	3200-1312	AO Drill Bit Ø2.9*250 with Block	1
13	3200-1313	AO Drill Bit Ø2.9*250 with Block	1
14	3200-1314	Drill & Tap Sleeve combined 2.5/3.5	1
15	3200-1315	Guider 2.0/4.0	1
16	3200-1316	Screwholding Forcep	1
17	3200-1317	Extraction Screw	1
18	3200-1318	Threaded Guider Wire Ø1.5*200mm	1
19	3200-1319	Threaded Guider Wire Ø1.5*200mm	1
20	3200-1320	K Wire Ø1.5*200mm	1
21	3200-1321	K Wire Ø1.5*200mm	1
22	3200-1322	Periosteal Dissector 12mm	1
23	3200-1323	Retractor 8.5mm	1
24	3200-1324	Retractor 15.5mm	1
25	3200-1325	Sleeve Key	1
26	3200-1326	Hollow	1
27	3200-1327	Periosteal Dissector 9mm	1
28	3200-1328	T-handle Quick Coupling	1
29	3200-1329	Tap for Cortical Screw 3.5mm	1
30	3200-1330	Tap for Cancellous Screw 4.0mm	1
31	3200-1331	Countersink	1
32	3200-1332	Tap for Locking Screw 3.5mm	1
33	3200-1333	Star Screwdriver T15	1
34	3200-1334	Star Screwdriver T15	1
35	3200-1335	Depth Gauge	1
36	3200-1336	Straight Handle Quick Coupling	1
37	3200-1337	Sharp Reduction Forcep (190mm)	1
38	3200-1338	Torque Handle 1.5N.m	1
39	3200-1339	Self-centering Bone Holding Forcep (190mm)	2
40	3200-1340	Oblique Reduction Forcep (170mm)	1
41	3200-1341	Bending Iron	1
42	3200-1342	Aluminium Box	1

Rib Reconstruction Locking Plate 3200-09





NO.	REF	Description	Qty.
1	3200-0901	Locking Sleeve	1
2	3200-0902	Locking Sleeve	1
3	3200-0903	Screwdriver	1
4	3200-0904	Screwdriver	1
5	3200-0905	Limited Drill Bit 2.2*8mm	1
6	3200-0906	Limited Drill Bit 2.2*10mm	1
7	3200-0907	Limited Drill Bit 2.2*12mm	1
8	3200-0908	Limited Drill Bit 2.2*14mm	1
9	3200-0909	Limited Drill Bit 2.2*16mm	1
10	3200-0910	Limited Drill Bit 2.2*18mm	1
11	3200-0911	Limited Drill Bit 2.2*20mm	1
12	3200-0912	Depth Gauge	1
13	3200-0913	Straight Handle	1
14	3200-0914	Straight Handle	1
15	3200-0915	Plate Holding Forcep	1
16	3200-0916	Plate Holding Forcep	1
17	3200-0917	Reduction Forcep	1
18	3200-0918	Reduction Forcep	1
19	3200-0919	Plate Bender	1
20	3200-0920	Plate Bender L	1
21	3200-0921	Plate Bender R	1
22	3200-0922	Drill Guider	1
23	3200-0923	Aluminium Box	1

Large Fragment Locking Plate Instrument Set (AO) 3200-14




NO.	REF	Description	Qty.
1	3200-1401	Hex Key	1
2	3200-1402	T-handle Quick Coupling	1
3	3200-1403	Transition Interface	1
4	3200-1404	Sleeve Key	1
5	3200-1405	Screwholding Forcep	1
6	3200-1406	Extraction Screw	1
7	3200-1407	Star Screwdriver T20	2
8	3200-1408	Tap for Cortical Screw 4.5mm	1
9	3200-1409	AO Drill Bit Ø3.2	2
10	3200-1410	Countersink	1
11	3200-1411	Wire Sleeve 2.0*5.0mm	1
12	3200-1412	Drill & Tap Sleeve combined 3.2/4.5mm	1
13	3200-1413	Locking Sleeve 5.0mm	1
14	3200-1414	Locking Sleeve 4.2mm for Thin Plates	2
15	3200-1415	Locking Sleeve 4.2mm for Thick Plates	2
16	3200-1416	Locking Sleeve 4.2mm for Thick Plates	1
17	3200-1417	Wire Sleeve 2.0*4.2mm	1
18	3200-1418	Wire Sleeve 2.0*4.2mm	1
19	3200-1419	Tap for Locking Screw 5.0mm	1
20	3200-1420	AO Drill Bit Ø4.2	1
21	3200-1421	AO Drill Bit Ø4.2 with Limitator Block	1
22	3200-1422	Cannulated Drill Bit for 6.5 Cannulated Locking Screw	1
23	3200-1423	K Wire Ø2.0*250mm	2
24	3200-1424	Threaded Guide Wire Ø2.0*250mm	2
25	3200-1425	Threaded Guide Wire Ø2.5*300mm	2
26	3200-1426	Periosteal Dissector 8mm	1
27	3200-1427	Sharp Reduction Forcep (200mm)	1
28	3200-1428	Cannulated Screwdriver for 6.5 Cannulated Locking Screw	1
29	3200-1429	Periosteal Dissector 15mm	1
30	3200-1430	Depth Gauge	1
31	3200-1431	Hollow	1
32	3200-1432	Bending Iron	1
33	3200-1433	Oblique Reduction Forcep (230mm)	1
34	3200-1434	Torque Handle 4.0N.m	1
35	3200-1435	Dissector 18mm/42mm	2
36	3200-1436	Self-centering Bone Hold Forcep(270mm)	1
37	3200-1437	Straight Handle Quick Coupling	1
38	3200-1438	Aluminium Box	1

Osteotomy Locking Plate Instrument Set 3200-08





NO.	REF	Description	Qty.
1	3200-0801	Quick Coupling Drill Bit Ø3.5*150	1
2	3200-0802	Quick Coupling Drill Bit Ø3.5*150	1
3	3200-0803	Quick Coupling Tap HC.0*220	1
4	3200-0804	Quick Coupling Drill Bit Ø4.3*280	2
5	3200-0805	Osteotome 10mm	1
6	3200-0806	Osteotome 15mm	1
7	3200-0807	Osteotome 20mm	1
8	3200-0808	Osteotome 25mm	1
9	3200-0809	Pointed Kirschner Wire Ø2.0*280	2
10	3200-0810	Pointed Kirschner Wire Ø2.5*280	2
11	3200-0811	Drill Bit KirschnerWire Ø2.5*300	2
12	3200-0812	Depth Gauge 0-120mm	1
13	3200-0813	Section Height Monitor	1
14	3200-0814	Angle Measurer	1
15	3200-0815	Torque Handle 4.0N.m	1
16	3200-0816	Spring Drill Ø3.5/4.3	1
17	3200-0817	Straight Quick Coupling Handle	1
18	3200-0818	Measuring Ruler	1
19	3200-0819	Quick coupling Hex Screwdriver SW3.5*100	1
20	3200-0820	Tap Guide Drill Sleeve Ø5.0*100	1
21	3200-0821	Guide Pin Guide Ø2.0	1
22	3200-0822	Guide Pin Guide Ø2.0	1
23	3200-0823	Drill-oriented Drill Sleeve Ø4.3*150	1
24	3200-0824	Drill-oriented Drill Sleeve Ø4.3*150	1
25	3200-0825	Aluminium Box	1

2.4/2.7 Multi-axial Locking Plates Instrument Set 3200-05



NO.	REF	Description	Qty.
1	3200-0501	Screw Rack	1
2	3200-0502	Multi-axial Cone Guide Ø2.1	1
3	3200-0503	Multi-axial Cone Guide Ø2.0	1
4	3200-0504	Depth Gauge	1
5	3200-0505	Torx Head Screwdriver T8	1
6	3200-0506	Multi-axial Straight Guide Ø2.1	1
7	3200-0507	Multi-axial Straight Guide Ø2.0	1
8	3200-0508	Screwdriver Handle	1
9	3200-0509	Drill Sleeve Ø2.0	1
10	3200-0510	Drill Sleeve Ø2.0	1
11	3200-0511	Drill Sleeve Ø2.1	1
12	3200-0512	Drill Sleeve Ø2.1	1
13	3200-0513	Plate Bender 2.4/2.7mm	1
14	3200-0514	Plate Bender 2.4/2.7mm	1
15	3200-0515	Drill Guide Ø2.1/2.7	1
16	3200-0516	Drill Guide Ø2.0/2.4	1
17	3200-0517	AO Drill Bit Ø2.0mm	1
18	3200-0518	AO Drill Bit Ø2.0mm	1
19	3200-0519	AO Drill Bit Ø2.1mm	1
20	3200-0520	AO Drill Bit Ø2.1mm	1
21	3200-0521	Drill Guide Ø2.1/2.7	1
22	3200-0522	Drill Guide Ø2.0/2.4	1
23	3200-0523	Torque Wrench 0.8N.m	1
24	3200-0524	Screw Holding Sleeve	1
25	3200-0525	Countersink Drill Ø2.4	1
26	3200-0526	Countersink Drill Ø2.7	1
27	3200-0527	Torx Head Screwdriver T8 (Quick-coupling)	1
28	3200-0528	Quick Coupling Tap HC2.4	1
29	3200-0529	Quick Coupling Tap HA2.4	1
30	3200-0530	Quick Coupling Tap HA2.7	1
31	3200-0531	Quick Coupling Tap HC2.7	1
32	3200-0532	Guide Pin	1
33	3200-0533	Guide Pin	1
34	3200-0534	Aluminium Box	1

FNS Instrument Set 3300-04




NO.	REF	Description	Qty.
1	3300-0401	Torque Limiting Attachment, 4 Nm, for AO Reaming Coupler	1
2	3300-0402	Large Cannulated Handle, Quick Coupling, 12 mm Hexagonal	1
3	3300-0403	Transition Interface	1
4	3300-0404	Protection Sleeve	1
5	3300-0405	SW3.5 Screwdriver Shaft for Locking Screw	1
6	3300-0406	T25 StarDrive Screwdriver Shaft for Bolts+Antirotation Screws	1
7	3300-0407	Screwdriver Shaft	1
8	3300-0408	Multifunction Rod, for Insertion Instruments	1
9	3300-0409	Depth Gauge up to 100 mm	1
10	3300-0410	Direct Measuring Device, for 3.2 mm Guide Wires	1
11	3300-0411	Insertion Handle	1
12	3300-0412	130 Degree Angled Guide, for 3.2 mm Guide Wires	1
13	3300-0413	Insert for Insertion Handle	1
14	3300-0414	Cylinder, for Insertion Instruments	1
15	3300-0415	Correction Guide, for 3.2 mm Guide Wires	1
16	3300-0416	Complete Opening Drill Bit/Reamer	1
17	3300-0417	3.2 mm Guide Wire, 400 mm	4
18	3300-0418	4.3 mm Drill Bi with stopper, 413 mm in Length	1
19	3300-0419	Aluminium Box	1

Sternum Locking Plate Instrument Set 3200-17



NO.	REF	Description	Qty.
1	3200-1701	Screw Candy	1
2	3200-1702	Temporary holder HC3.0	1
3	3200-1703	Guide drill 2.3	2
4	3200-1704	Temporary holder HC3.0	2
5	3200-1705	Quadrangle screwdriver SW1.7	1
6	3200-1706	2.3 Limited Drill Bit 6mm	1
7	3200-1707	2.3 Limited Drill Bit 8mm	1
8	3200-1708	2.3 Limited Drill Bit 10mm	1
9	3200-1709	2.3 Limited Drill Bit 12mm	1
10	3200-1710	2.3 Limited Drill Bit 14mm	1
11	3200-1711	2.3 Limited Drill Bit 16mm	1
1.2	3200-1712	Plate bender L=200mm	2
13	3200-1713	Plate shears L=190mm	2
14	3200-1714	Plate Holding Compression Spreader HC3.0	1
15	3200-1715	Percutaneous Guiding Drill 3.0	1
16	3200-1716	Depth Gauge 0-30mm	1
17	3200-1717	Straight Handle	1
18	3200-1718	Aluminium Box	1

Subtalar Screw Instrument Set 430-004



NO.	REF	Description	Qty.
1	430-760	Peeler	1
2	430-800	Cannulated Screwdriver SW4.0	1
3	430-801	Cannulated Screwdriver SW4.5	1
4	430-740	Trail Φ 5.5 $\times$ 21	1
5	430-741	Trail Φ 6.5 $\times$ 21	1
6	430-742	Trail Φ 7.3 $\times$ 21	1
7	430-743	Trail Φ 8.5 $\times$ 21	1
8	430-744	Trail Φ 9.5 $\times$ 21	1
9	430-750	Trail Φ 8.5 $\times$ 15	1
10	430-751	Trail Φ 10 $\times$ 15	1
11	430-752	Trail Φ 11 $\times$ 15	1
1.2	430-753	Trail Φ 12 $\times$ 15	1
13	411-730	Bone Hook	1
14	411-730	Bone Hook	1
15	430-780	Screw Holding Forcep	1
16	430-790	Bone Cutter	1
17	101-247	Rechat Handle	1
18	430-770	K-wire 2mm	2
19	430-810	Alumium Box	1

Small Fragment Locking Screw Rack 1000-0107



Size:40*14*11cm

3.5mm Locking Screw Length:12-95mm

3.5mm Cortex Screw Length:12-60mm

4.0mm Cancellous Screw Length:12-50mm

Large Fragment Locking Screw Rack 1000-0108




Size:35*13*12cm

5.0mm Locking Screw Length:22-95mm

4.5mm Cortex Screw Length:22-60mm

6.5mm Cannulated Locking ScrewLength:50-110mm

Radius Plate Box 1000-0128



For 3-7 holes VA Radius Locking Plate and
Locking Screw